GeneDisc® Rapid Microbiology System

A revolutionary step in advanced microbial quality monitoring

Confidential

ISGA October 2011





Disclaimer

This presentation is the Confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer.





A multi-parametric food pathogen testing system

- Easy to use, no need for strong technical skills
- Flexible, modular system, scale to your test volume
- Real time PCR
- At-a-glance display of results, enables rapid decisions
- LIMS connectivity for data management
- Salmonella and E.coli test results on the same sample



Real time PCR Overview

- All microorganisms contain unique DNA sequences (targets)
- After DNA Extraction,





Outline Description of the System

- There are three core components to the system:
 - The Extraction Pack and associated consumables
 Universal method for microbial DNA extraction in all food matrices

Confidential

- The GeneDisc plate and associated consumables
- To assay extracted microbial DNA
- The GeneDisc Cycler

To run assays to detect or quantify this DNA







Cycler: the qPCR assay instrument

GeneDisc[®] Real time PCR technology



ife Sciences

- Assay Consumable
 - GeneDisc[®] Plate
 - TaqMan technology
 - Primer/Probes preloaded
 - Inhibition control/sample
 - Negative control/sample
 - Plate for 6 or 12 samples
 - Bar code for traceability
 - 1 tube of master mix / plate
 - Minimal volume to work with 20 μL
 - Limited risk of cross contamination

GeneDisc® Real time PCR technology



Life Sciences

PAL

- Real-Time PCR instrument
 - GeneDisc[®] Cyclers

Designed for routine analysis

- Adaptable, from 1 to 8 units
- 1 to 8 GeneDisc plates simultaneously, ie. 1 to 96 samples



Process flow

Enrichment & DNA extraction

Pall Food DNA Extraction kit

Screening by real time PCR

GeneDisc & GeneDisc Cycler



GeneDisc Cycler Heating Process :



PCR Operations with GeneDisc Plate





Barcode Reading

- GeneDisc bar-coded;
- Barcode selects appropriate qPCR parameters.

Genel



GeneDisc Filling



GeneDisc Loading – PCR • run starts.

- 1 Master Mix tube ↔ 1 GeneDisc -Limited risk of contamination;
- GeneDisc filled (vacuum system), placed into GD Cycler, lid closed
- Assay starts automatically no further operator intervention
 - Disc completely sealed, minimising contamination risk.

PALL Pall Corporation	GeneDisc [®] Cycler			1.0.1.02 st be (TECH)	
	Analysis report on (Cycler 2			
Cycler 1	User	st bo	Lot		010010111321ST
Analysis	Date	ate 04/01/11 03:10 PM W		rning(s)	2 alarm(s) 🔻
	PCR application	GeneDisc® Shi	gaToxic E. coli	06_1	
	General By sample	Curves By Well	s Notes V	Validation	
Records Cycler 3	Sector	Sample		Interpretatio	n
	1 s1		Presence of pathoger	nic E. coli O157	
Cycler 4	2 82		Absence of pathogenic STEC		
	4 54		Absence of pathogen	ic STEC	
Cyclet 5	5 s5		Presence of pathoger	nic STEC	
	6 56		Presence of pathoger	nic STEC	
2					
Help					
Cycler 7					
Log off Cycler 8					
				\checkmark	
		Confidential 11			







PALL



ALL Pall Corporation	GeneDisc [®] Cycler			1.0.0.9 demo demo (ADMIN)	
	Analysis report				
	Parameter Preview				
Analysis	Template 1 🔻	Sample	5:		
	Analysis lab. address				
	Client reference	La dela constante La dela const			
	Device serial number	Brit Bud	Over advess		
Records	GeneDisc lot No.	<u>Citet reference</u> Analysis performed the 1977/1984	(2) Pe		
	Master Mix lot No.	PD estudios - Sendinal for Sergie	dragens & sol (CP B.)		
	Notes	Matter Wester Ro. 1			
Admin.	Curves	Target mult			
	Øct	L Gr Kan In , Sin In Con			
2		<u>Saturbanan</u> Prosence: dinan patricipani, 6, a	u-Şeşî		
Help			Agentical Tay I Same James Al Braz The 2019 19 (1922) All		
2			page 1		
Log off 10:22 AM			٢	0	
Life Sciences	Con	fidential			

PALL Pall Corporation	GeneDisc [®] Cycler			
	Analysis report			
2	Parameter Preview			
Analysis	Template 1	Sample	5:	
	Analysis lab. address	(Machine La Statement		
- E	Client reference Device serial cumber	Fair Constitutions 1 Gas do Consti UK/20 Binar		
Records	GeneDisc lot No.	Outstance	Olevit address	
	Master Mix lot No.	Analysis performant inst, 1011-101 POR application - Demokrati P	na an a	
	Notes	Manu Mular Au		
Admin.	Curves	199	- Nagathar control	
	₽ ct	- 190 - 190 - 190 - 190		
2		138		
lielp		1 100		
			Algermanility during during Al Brazi Ting (20 + 6.10) 10.000 ABB Digital and	
			Control 2	
Log off				
10:22 AM			\bigcirc	00
Life Sciences	Cor	nfidẹntial		

STEC: The example of O157



Design of the GeneDisc STEC & E. coli O157

- GeneDisc
 - Simultaneously detection of *E. coli* O157, STEC & Salmonella spp

Primers and probes design used for pathogenic genes are recommended by CEN ad hoc group for their high specificity

		Well #	FAM detection	ROX detection
		1	Negative CTRL	Inhibition CTRL
	2	2	iroB (Salmonella)	stx1/stx2*
	3	3	rfbE ₀₁₅₇ *	eae*
		*Primers a "STEC ad Standardi	and probes used are t hoc group" of the Eu zation (CEN) for a fut	hose proposed by the ropean Committee for ure ISO standard
PALL Life Sciences	*Perelle Confidential *Nielsen &	S., et. al., Andersen	Molec. Cell. Probes , J. Clin. Microbiol.,	3 18: 185-92 (2004) 41:2884-93 (2003)

Focus on Pathogenic *E.coli* O157 and non-O157 STEC



- The additional primers and probes (O121 & O45) are those recommended by USDA ARS (Dr Pina Fratamico)
- Added O104:H4 plate for testing related to outbreak in Germany







Chronology of E.Coli O104:H4 project

- 1. End of May 2011 \rightarrow <u>Outbreak started</u> being reported to media
 - Important pressure on German and EU authorities;
 - No official analytical method available for the outbreak strain, E. coli O104:H4.
- 2. Former study run with GD technology had some available sequences
 - Work closely to the <u>French Food Safety Agency (Anses)</u> for developing the right targets;
 - Adaptation on the GD technology and production of the first batch for quick validation by Anses.
- 3. By June 2011, German lab (Bfr Berlin managed by Dr Lothar Beutin) got the GD 0104:H4 for food sample screening
 - Investigations of 100's of samples;
 - Food survey identified <u>fenugreek seeds as the most probable source of infection</u>.
- 4. By end of June 2011, Pall was used for the <u>French outbreak investigations linked to *E. coli* O104:H4.</u>



Pall GeneDisc Technologies – STEC O104:H4



One Step → detection of various virulence gene and O104:H4 genes

- stx2
- wzx O104
- flic-H4
- *ter* (tellurite resistance)
- aggR (transcriptional regulator)
- Inhibition and negative control per sample

Filtering Sprout Irrigation Water

- Test more water without running more tests
- Increase sample size from 1 liter to 5 or 10 liters by using a concentration step. Direct flow – enrich filters





Filtering Sprout Irrigation Water

- Test more water without running more tests
- Increase sample size from 1 liter to 5 or 10 liters by using a concentration step. Tangential Flow – enrich concentrate.



Irrigation Water Challenges

Particle Size Distribution



Jonathan's Sprout Irrigation Water



Irrigation Water Challenges

Irrigation Water Concentration Factors				
Factor	Direct Flow	Tangential Flow		
Time	Sizing	Sizing		
Cost	Lower	Higher		
Cross Contamination	Single Use	Sanitize and Reuse		
Testing	Material Retained on Filters	Concentrate		



Thank you

Any Questions ?

Contact: Wayne Miller Email: <u>Wayne Miller@pall.com</u> Phone: 774-276-5854

An innovative step in advanced microbial quality monitoring

_ife Sciences

