



Pall GeneDisc[®] Technologies

From regulation to diagnostics

Simplified compliance to the new
European Regulation for sprouted
seeds

This presentation is the Confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer

© 2013 Pall Corporation

Eric Samuels

ISGA convention, april 22nd, 2015

Rotterdam, the Netherlands



From Regulation to Diagnostocs

Content:

- Intro Pall GeneDisc Technologies
- EU-Regulation: today and tomorrow
- GeneDisc Solution for STEC screening



History

- French Biotech Company, GeneSystems, founded in 2000.
- Product Offering: Molecular diagnostic solution provider in the field of food analysis, environment and forensic science.
- Proprietary technology in molecular biology (PCR instrument, PCR reagents, DNA extraction and purification).
- Acquired by Pall Corporation in September 2008. Product development and placement across multiple Pall divisions: Medical and Food & Beverage.





A European regulation for sprouts !

And a new guidance for all foods?



European Food Safety Regulation – Context

Food Business Operators (FBO) are concerned, among others, by the following EC regulations:

- EC regulation 852/2004 which sets FBO general rules on hygiene of foodstuffs (focus on HACCP)
- EC regulation 2073/2005 which sets harmonized criteria for certain micro-organisms and how to perform tests (e.g. *Salmonella* criteria is listed for sprouted seeds)

<http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2005:338:0001:0026:EN:PDF>



European Food Safety Regulation – Reasons for evolution

Several factors have occurred recently and raised concerns about Shiga-Toxic *Escherichia coli* (STEC) in sprouts:

- *E. coli* (STEC) outbreak in EU (May 2011; 3916 Human cases and 47 deaths)
- EFSA Scientific Opinion on the risk posed by STEC in seeds and sprouted seeds

Conclusion: Need to extend microbiological criteria to STEC for the sprouted seeds production chain



European Food Safety Regulation – Evolution

Following these events, Commission Regulation 209/2013 Amending EC regulation 2073/2005 has been published

12.3.2013

EN

Official Journal of the European Union

L 68/19

COMMISSION REGULATION (EU) No 209/2013

of 11 March 2013

amending Regulation (EC) No 2073/2005 as regards microbiological criteria for sprouts and the sampling rules for poultry carcasses and fresh poultry meat

(Text with EEA relevance)

<http://eur-lex.europa.eu/JOHtml.do?uri=OJ:L:2013:068:SOM:EN:HTML>



What is Commission Regulation 209/2013?

This Commission Regulation establishes:

- A new microbiological testing requirement for sprouts which concerns STEC (serogroup O157, O26, O103, O111, O145 and O104:H4)
- A method for this testing (ISO TS 13136)
- An acceptance criteria (Absence in 25g)
- Sampling rules:
 - Preliminary testing on a representative sample of all batches of seed
 - Testing at the stage where probability of finding STEC and *Salmonella* is the highest
 - Alternative testing of spent irrigation water is allowed
 - Frequency must be at least once a month

FBO can only place sprouts on the market if the results comply with regulation acceptance criteria for *Salmonella* and STEC



Sprouts testing regulation after amendment

Food Category	Micro-organisms	Sampling plan	Limits	Analytical reference method	Stage where the criterion applies
1.18 Sprouted seeds (ready-to-eat)	<i>Salmonella</i>	n=5; c=0	Absence in 25g	EN/ISO 6579	Products placed on the market during their shelf-life
1.29 Sprouts	Shiga toxin producing <i>E.coli</i> (STEC) O157, O26, O111, O103, O145 and O104:H4	All new batches n= 5; c= 0	Absence in 25g	EN/ISO 13136	Products placed on the market during their shelf life



And a new guidance for all foods?



EUROPEAN COMMISSION
HEALTH AND CONSUMERS DIRECTORATE-GENERAL

DRAFT

GUIDANCE DOCUMENT ON THE APPLICATION OF ARTICLE 14 OF REGULATION (EC) N°178/2002 AS REGARDS FOOD WHERE SHIGA TOXIN- PRODUCING ESCHERICHIA COLI (STEC) HAS BEEN DETECTED

- Almost ready for publication after acceptance by the EC


This document has been established for information purposes only. It has not been adopted or in any way approved by the European Commission.

The European Commission does not guarantee the accuracy of the information provided, nor does it accept responsibility for any use made thereof. Users should therefore take all necessary precautions before using this information, which they use entirely at their own risk.

- In the Netherlands implemented per 22 may 2014




A new guidance for all foods?





A new guidance for all foods?





What is ISO TS 13136?

The screenshot shows the ISO website's search interface. At the top, there's a navigation bar with links: Standards, About us, Standards Development, News, and Store. Below this, a search bar contains the text 'STEC'. To the right of the search bar are filters for 'Standards (1)', 'Publications & e-products (0)', and 'Site (1)'. A 'Search' button is next to the search bar. Below the search bar, there's a 'Sort by:' dropdown menu set to 'Relevance', and a link to 'Advanced Search'. Further down, there's a section titled 'Standards to be displayed:' with three checkboxes: 'Published' (checked), 'Under development' (unchecked), and 'Withdrawn' (unchecked). Below this, there's a link to 'Refine your query' and a result count '1 - 1 of 1'. The search result is for 'ISO/TS 13136:2012', which is marked with a checkmark. The description of the standard is: 'Microbiology of food and animal feed -- Real-time polymerase chain reaction (PCR)-based method for the detection of food-borne pathogens -- Horizontal method for the detection of Shiga toxin-producing Escherichia coli (STEC) and the determination of O157, O111, O26, O103 and O145 serogroups'. A link to 'More details »' is provided at the bottom of the result.

This technical specification published in November 2012 is relative to STEC testing and is based on Real Time PCR principle.

It describes :


- Enrichment protocol,
- DNA extraction,
- Screening for virulence factors (*stx* and *eae*),
- Identification of serogroups (O157, O111, O26, O103 and O145),
- Isolation.

O104:H4 is not included in the ISO TS 13136, but testing should be performed according to this 5 steps method for this serotype.




Flow-diagrams of ISO TS 13136 operating workflow

Flow-diagram of the screening procedure



Flow diagram of the isolation procedure





GeneDisc Solution



GeneDisc Overview

- Pall's GeneDisc system is a robust platform which allows accurate detection of microorganisms using the Real-Time Polymerase Chain Reaction (RT PCR) method
- The principle of Real-Time Polymerase Chain Reaction is the amplification of a specific DNA sequence





GeneDisc Cyclers

The GeneDisc Cycler is an high throughput real time PCR instrument for the detection of microorganisms.

GeneDisc Cycler features:

- High Throughput - test up to 96 samples each hour
- Flexible and Modular - up to eight individual Real-Time PCR units can be used at once
- Rapid - Run time less than 1 hour
- Easily fits in lab - Small footprint
- Fast decision making - Real time measurement with at a glance results
- Seamless data transfer - LIMS compatible





GeneDisc Plate

GeneDisc Plate features:


- Single/simultaneous detection of different organisms
- Samples to test results in hours
- Easy-to-use system, minimal operator intervention required
- Barcode control of consumables and assay parameters
- Key reagents pre-loaded in a unique, sealed plate (no operator handling of primers...)
- Unique design allows automatic filling of the plate
- High specificity, reliability and reproducibility of run





GeneDisc overview

Enrichment and DNA Extraction From sample to DNA



GeneDisc Cyclor From DNA to result



GeneDisc Plate Ready to use PCR plate




GeneDisc Method for STEC

- GeneDisc method follows ISO TS 13136 requirements in term of target gene detection and workflow :

1. Screening for virulence factors + O-groups

- Use *E.coli* TOP7 GeneDisc Plate



GTOP7EC106006


Well	FAM Detection	ROX Detection
1	O145	Inhibition control
2	Vir. O45-O103-O121	stx1-2
3	Vir. O26	Vir. O145-O157
4	Vir. O111	O111
5	O157	O26
6	-	O45-O103-O121

GeneDisc Method for STEC

- GeneDisc method follows ISO TS 13136 requirements in term of target gene detection and workflow :

1. Screening for virulence factors (* eu guidance)

- Use ShigaToxic *E.coli* GeneDisc Plate



1
2
3


Well Number	FAM Fluorescence Dye Detection	ROX Fluorescence Dye Detection
1	PCR Negative control	PCR inhibition control
2	<i>stx2</i>	<i>stx1</i>
3	<i>E. coli</i> O157	<i>eae</i>



GeneDisc Method for STEC

2. Identification of Top 7 (O26, O103, O111, O145) if screening for virulence factors is positive

➤ Use EHEC 5 ID GeneDisc Plate



Well Number	FAM Fluorescence Dye Detection	ROX Fluorescence Dye Detection
1	PCR Negative control	PCR Inhibition control
2	<i>E. coli</i> O145	
3	Flagellar H7	<i>E. coli</i> O111
4	Flagellar H7	<i>E. coli</i> O111
5	<i>E. coli</i> O26	<i>E. coli</i> O103
6	<i>E. coli</i> O26	<i>E. coli</i> O103


3. Screening for O104:H4 when *stx2* is positive

➤ Use *E. coli* O104:H4 GeneDisc Plate




ISO TS 13136 and GeneDisc Solution

Flow-diagram of the screening procedure



Flow diagram of the isolation procedure





A solution to your needs



GeneDisc STEC solution

- Meet EN 2073 – All required testings performed at once :
 - Detection of all pathogenic STEC Top 7 and O104:H4
 - Simultaneous detection of STEC and *Salmonella* spp. available
- The method has been tested on seeds, sprouts and water.



GeneDisc STEC solution

- This solution enables an easy and accurate testing for STEC
- This solution is compliant to ISO/TS 13136 requirements
- This solution minimal handling of samples allows to easily perform STEC testing.
- GeneDisc method for STEC screening and serogroup identification has been AOAC approved





A solution to your needs

Your needs	Our solution
✓ Brand protection/Compliance to Regulation	✓ STEC testing compliant to ISO TS 13136 specifications ✓ <i>Salmonella</i> testing validated alternative method to ISO 6579
✓ Minimal handling	✓ Simultaneous STEC/ <i>Salmonella</i> detection (one enrichment, one sample preparation, one run)
✓ Earlier batch release	✓ Rapid – Result obtained in as fast as 20h
✓ Adaptive to your throughput and sampling plan	✓ Modular system
✓ Smooth implementation	✓ Ease of use



Further information

Several external accredited Laboratories are already equipped with the the GeneDisc Solution for STEC screening.

For further information, please contact us at:

genedisc@pall.com

Or:

Eric_Samuels@europe.pall.com

+31 20706121

Please, visit us at:

www.pall.com/genedisc



Dear Customer,

Recently the ISO/TS 13136:2011 has been published:

ISO/TS 13136:2011 (E)

Microbiology of food and animal feed — Real-time polymerase chain reaction (PCR)-based method for the detection of food-borne pathogens — Horizontal method for the detection of Shiga toxin-producing Escherichia coli (STEC) and the determination of O157, O111, O26, O103 and O145 serogroups

Microbiologie des aliments — Méthode basée sur la réaction de polymérisation en chaîne (PCR) en temps réel pour la détection des micro-organismes pathogènes dans les aliments — Méthode horizontale pour la détection des Escherichia coli producteurs de Shigatoxines (STEC) appartenant aux sérogroupes O157, O111, O26, O103 et O145

Pall GeneDisc Technologies SA certify that the primers and probes used in GeneDisc products listed in attachment A, are as proposed in the ISO TS 13136 for the stx1, stx2, eae, O157, O111, O103, O145, O26.

For the target genes that are patented, Pall GeneDisc Technologies SA has obtained a licence agreement from ANSES.


Therefore, GeneDisc PCR Solution is fully compliant with the primer and probes sets described in the ISO/TS 13136.


Further our GeneDisc Solutions for STEC screening has been previously approved by AOAC-PTM against the first draft of the ISO/TS 13136.

We kindly suggest you add this information to your file and hope to work closely with you on the detection of STEC in food products.

With kind regards,

Pall GeneDisc Technologies SA


Mrs Sirine Assaf
Product Manager GeneDisc


Mrs Nathalie Kerriguy
QA manager GeneDisc



www.pall.com/green

This presentation is the Confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer

© 2013 Pall Corporation